



# Temporary Pacemakers



Medtronic, Inc., Minneapolis, MN  
June 2007

# Model 5388 Dual Chamber Temporary Pacemaker

1. Pace/Sense LEDs
2. Lock/Unlock Key
3. Lock Indicators
4. Rate Dial
5. Atrial Output Dial
6. Ventricular Output Dial
7. Menu Parameter Dial
8. Parameter Selection Key
9. Menu Selection Key
10. Pause Key
11. Power On Key
12. Power Off Key
13. Emergency/Asynchronous Pacing Key
14. Lower Screen
15. Ventricular Output Graphics
16. Atrial Output Graphics
17. Upper Screen
18. Rate Graphics
19. Setup Indicators
20. DDI Indicator
21. Low Battery Indicator
22. Setup Labels




Medtronic, Inc., Minneapolis, MN  
June 2007

# Off / On Keys

## Values at Power-On

### Dual Chamber Pace/Sense

- **RATE** 80 ppm
- **UPPER RATE** 110 ppm



**OFF**

Push twice


**ON**

Push once



Medtronic, Inc., Minneapolis, MN  
June 2007

# Emergency Key




Always available – Single key press enters Emergency mode

## Emergency Pacing Values

- |               |              |
|---------------|--------------|
| • RATE        | Current Rate |
| • A OUTPUT    | MAX          |
| • V OUTPUT    | MAX          |
| • PACING      | ASYNC        |
| • NO SENSING! |              |


Use caution when setting the device to asynchronous modes.

# Pause Key – Check Patient's Intrinsic Rhythm



# Indicators

Pace/Sense Indicators  
(Flashing Lights)




Pace/Sense  
Setup Indicators  
("how is the device setup?")

Padlock  
Indicator  
("is the device  
locked?")



# Lock / Unlock



Locks Rate, V Output, A Output dials

Lock/Unlock Key

Lock Indicator

Padlock  
Icon


Flashing Key  
Icon

**To unlock push the  
“Lock/Unlock” Key**

**Emergency key is  
always available**

# Rate and Output Adjustments

## *Single or Dual Chamber Pacing With Only 3 Dials!*



Rate Dial

Max rate of  
**200bpm** for  
pediatrics

Atrial Output Dial

For Single Chamber pacing, turn  
OFF Atrial output

Ventricular Output Dial



# Lower Screen Menus

**DDD**

A Sensitivity **0.5 mV**

V Sensitivity **2.0 mV**  
 20 10 0.8

A-V Interval **170 mS**

A Tracking **ON**

1

Menu 1: Pacing Parameters

**DDD**

Upper Rate **110 ppm**  
 80 155 230

PVARP **300 mS**

A-V Interval **170 mS**

SETTING

2

Menu 2: Rate-Based Pacing Parameters

**DDD**

RAP **320 ppm**  
 80 440 800

Press **SELECT**

to DELIVER

RAPID ATRIAL PACING

3

Menu 3: Rapid Atrial Pacing

**AAI**

DIAL-A-MODE **DOO**  
 DDD.....DVI.....DOO.....VVI

Press **SELECT**

to ACTIVATE


new Mode

M


Menu M: Dial-A-Mode

# Device Usage - Cable Connectors


- Connector pins on the lead(s) must be fully inserted in the patient connector block
- Observe polarity
- **Finger tighten only** – no tools!



# Cable to Device Connections



# Emergency Connections




Temporary use only -  
Leads do **not** lock in place

# Battery Operation

*Replace the battery for each new patient*

## Low Battery Indicator



### Low Battery Indicator

When indicator first appears you have approximately **24 hours** of pacing remaining

### Check Status



Check battery status at least **twice** per day – replace battery when indicator is on

Replace battery at least once per week when device is in continuous use



Medtronic, Inc., Minneapolis, MN  
June 2007

# Battery Replacement



**15 seconds** of pacing provided  
while changing 9V battery  
(note: battery polarity is reversible)



Medtronic, Inc., Minneapolis, MN  
June 2007



# Pacing Setup Table

**Model 5388 Pacing Setup Table**

	<b>AOO*</b>		<b>VOO</b>		<b>AAI</b>		<b>VVI</b>		<b>DOO</b>		<b>DVI</b>		<b>DDD</b>		<b>DDI</b>	
<b>Setup Indicators</b>	PAGE	SENSE	PAGE	SENSE	PAGE	SENSE	PAGE	SENSE	PAGE	SENSE	PAGE	SENSE	PAGE	SENSE	PAGE	SENSE
	A		V		A	A	V	V	A+V		A+V	V	A+V	A+V	A+V	A+V
<b>Instructions</b>																
<b>1. Set Output</b>																
<b>A Output</b>	On		Off		On		Off		On		On		On		On	
<b>V Output</b>	Off		On		Off		On		On		On		On		On	
<b>2. Set Sensitivity</b>																
<b>A Sensitivity</b>	ASYNC		NA		On		NA		ASYNC		ASYNC		On		On	
<b>V Sensitivity</b>	NA		ASYNC		NA		On		ASYNC		On		On		On	
<b>3. Set</b>																
<b>A Tracking</b>	NA		NA		NA		NA		NA		NA		On		Off	
<b>*Caution:</b> DAD and OOO are accessible modes, but are not recommended. Refer to "Controls, Indicators, and Other Features" in the technical manual.														NA: Not Applicable		

# NBG Codes

## 1st Letter

### Chamber(s) Paced

A = atrium

V = ventricle

D = dual (both atrium  
and ventricle)

## 2nd Letter

### Chamber(s) Sensed

A = atrium

V = ventricle

D = dual

O = none

## 3rd Letter


### Response to Sensing

I = inhibit  
(Demand mode)

T = triggered


D = dual

O = none (Asynch)



# Setting Sensitivity

The degree that the pacing system “sees” or senses signals, controlled by the sensitivity setting which is graduated in millivolts (mV)



Sensitivity settings should provide at **least a 2:1 safety margin**

# Sensing Threshold Procedure

1. Set rate at least 10 ppm below patient's intrinsic rate.
2. Adjust output: Set OUTPUT to 0.1 mA (A OUTPUT for atrial threshold; V OUTPUT for ventricular threshold).
3. Highlight SENSITIVITY (atrial or ventricular) (Menu 1).
4. Decrease SENSITIVITY: Slowly turn **MENU PARAMETER** dial counterclockwise until pace indicator flashes continuously.
5. Increase SENSITIVITY: Slowly turn **MENU PARAMETER** dial clockwise until sense indicator flashes and pace indicator stops flashing. *This value is the sensing threshold.*
6. Set SENSITIVITY to half (or less) the threshold value. *This provides at least a 2:1 safety margin.*
7. Restore RATE and OUTPUT to previous values.



# Atrial/Ventricular Sensing Thresholds

Sensing




Atrial Undersensing



# Capture

Depolarization of cardiac muscle  
following an electrical stimulus





# Stimulation Threshold

The **minimum** output pulse needed to consistently capture the heart



Set output to 2-3 times stimulation threshold

# Stimulation Threshold Procedure

1. Set RATE at least 10 ppm above patient's intrinsic rate.
2. Decrease OUTPUT: Slowly turn **OUTPUT** dial counterclockwise until ECG shows loss of capture.
3. Increase OUTPUT: Slowly turn **OUTPUT** dial clockwise until ECG shows consistent capture.  
*This value is the stimulation threshold.*
4. Set OUTPUT to a value **2 to 3 times greater** than the stimulation threshold value.  
*This provides at least a 2:1 safety margin.*
5. Restore RATE to previous value.




# Atrial/Ventricular Stimulation Thresholds

## Capture



## Loss of Ventricular Capture




# Model 5348 Single Chamber Pacemaker



Medtronic, Inc., Minneapolis, MN  
June 2007

# Basic Device Operation

## Turning device off



Push both buttons at the **same** time to turn OFF

## Battery Replacement

Push Buttons at the same time





# Locking Feature and Emergency

## Locking Feature

Slide plastic cover over dials to protect against changes in settings

## Emergency Pacing

- Rate                      Check Rate
  - Output                  Turn to **MAX**
  - Sensitivity              Turn to **ASYNC**
- Use caution when setting the sensitivity to asynchronous*





# Troubleshooting Pacemaker Performance



Medtronic, Inc., Minneapolis, MN  
June 2007


# Troubleshooting Process

1. Gather information
2. Identify the problem and possible cause
3. Identify the solution and carry out corrective procedures

# Loss of Capture

Electrical stimuli delivered by the pacemaker does not initiate depolarization of the atria or ventricle

# Loss of Capture



## Possible Causes

- Threshold rise
- Fractured/dislodged lead
- Battery depletion
- QRS not visible
- Tissue is refractory
- Faulty cable connections

## Corrective Measures

- Increase output (mA)/check thresholds
- Replace/reposition lead
- Replace battery
- Adjust ECG
- Assess mode selection
- Check connections
- Switch polarity (epicardial system)




Medtronic, Inc., Minneapolis, MN  
June 2007

# No Output

Pacemaker fails to emit stimuli  
at the programmed intervals

# No Output



## Possible Causes

- Battery depletion
- Pacemaker OFF
- Faulty cable connections
- Fractured/dislodged lead
- Oversensing

## Corrective Measures

- Replace battery
- Verify pacemaker settings
- Check cable connections
- Replace/reposition lead
- Verify/adjust sensitivity




Medtronic, Inc., Minneapolis, MN  
June 2007



# Undersensing

Failure of the pacemaker to  
sense  
intrinsic R-waves or intrinsic  
P-waves

# Undersensing



## Possible Causes

- Decreased QRS voltage
- Fractured/dislodged lead
- Battery depletion
- Inappropriate sensitivity setting
- Fusion beat

## Corrective Measures

- Increase sensitivity
- Replace/reposition Lead
- Replace Battery
- Sensing test/increase sensitivity



# Fusion/Pseudofusion Beats



Intrinsic Beat



Paced Beat



Fusion Beat



Pseudofusion Beat

# Oversensing

Inhibition of the pacemaker by  
events  
pacemaker should ignore, e.g. EMI,  
T-waves and myopotentials

# Oversensing



## Possible Causes

- Fractured/dislodged lead
- Environmental interference
- T-wave oversensing
- Faulty cable connections

## Corrective Measures

- Replace/reposition lead
- Eliminate interference
- Sensing test/decrease sensitivity
- Check connections

# References

- Dubin D. Rapid Interpretation of ECGS. 6<sup>th</sup> ed. Tampa: Cover, Inc; 2000.
- Ellenbogen KA, Wood MA. Cardiac Pacing & ICDS. 3<sup>rd</sup> ed. Malden: Blackwell Science, Inc; 2002.
- Fogoros RN. Electrophysiologic Testing. 3<sup>rd</sup> ed. Malden: Blackwell Science, Inc; 1999.
- Hayes DL, Lloyd MA, Friedman PA. Cardiac Pacing and Defibrillation: A Clinical Approach. Armonk: Futura; 2000.
- Moses HW, Moulton KP, Miller BD, et al. A Practical Guide to Cardiac Pacing. 4<sup>th</sup> ed. Boston: Little, Brown; 1995.



Medtronic, Inc., Minneapolis, MN  
June 2007